

Polar Equine Electrode Set Polar Equine T54H (Wearlink 31C) Polar Equine T56H (Wearlink W.I.N.D.)


1

2

3

4

ENGLISH

How to attach the Electrode Set, T54H (Wearlink 31C) or T56H (Wearlink W.I.N.D.).

First of all, identify the 2 electrode pads on each side of the transmitter. The electrode with velcro on the back is the positive electrode and logically the other electrode is the negative electrode. In the middle, you'll find the transmitter pocket.

The transmitter picks up very small electrical impulses emitted by the heart and for the transmitter to read the heart rate properly, we need to insure a good contact to the horse's skin. Make sure you wet the electrodes properly, the integrated sponge will keep up some moisture, but you still need to wet to horse's coat, especially as your horse is not yet sweating. If your horse has a thick winter coat, you can shave the spot where you will place the electrodes, this will greatly improve the heart rate signal quality.

1) Place the positive electrode on the left or right side of the withers under the saddle. Ensure that the electrode is flat against the horse's skin. The rider's weight will keep the electrode in place.

2) Attach the WearLink 31C or W.I.N.D. transmitter into the transmitter pocket and close it. This pocket protects your transmitter from shocks and scratches and prevents it from falling off.

3) Attach the transmitter to the saddle with one of the rubber strap.

4) Attach the negative electrode to the girth using a rubber strap and tighten the girth. If you are riding with a loose girth, you might want to place a sponge behind the negative electrode to insure contact with the skin.

Start your receiver (refer to the receiver user manual for instructions) and hold it next to the transmitter, within a maximum range of 80cm-1m (WearLink 31C) or 8-10m (WearLink W.I.N.D.).

The horse's heart rate will appear on the receiver within 5 to 10 seconds. Check that you are not near other people with heart rate monitors or any source of electromagnetic disturbances, like high voltage power lines, overhead lines of electric railways, electric bus lines or trams, car motors, some motor driven exercise equipment or cellular phones.

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