

POLAR®

POLAR®

Manufactured by
Polar Electro Oy
Professorintie 5
FIN-90440 KEMPELE
Tel +358 8 5202 100
Fax +358 8 5202 300

www.polar.fi

17923224.00 ENG/FRA/ESP/PRT A



USER MANUAL
MANUEL D'UTILISATION
MANUAL DEL USUARIO
MANUAL DO USUÁRIO

F92ti™ Heart Rate Monitor
Cardiofréquencemètre F92ti™
Monitor de Ritmo
Cardíaco F92ti™
Monitor de Freqüência
Cardíaca F92ti™

Congratulations on purchasing a Polar Heart Rate Monitor. We are proud to assist you in achieving your personal fitness goals. This guide will provide you with the basic information on how to exercise with your Polar Heart Rate Monitor and how to get the maximum benefit from its functions.

To ensure the continuous improvement of Polar products, please complete and mail the enclosed Customer Registration Card within two weeks of the date of purchase.

Thank you for choosing Polar.

Visual Cues to Read this Guide



Indicates important issues. Please, pay attention!



Indicates a useful tip. It is good to know these things.



Indicates a tip to do things faster.

Important to Remember!

- Adjust the strap length snug enough and wet the Polar transmitter electrodes to ensure flawless heart rate measurement.
- Do not operate the buttons of the wrist unit under water.
- Wash and dry your transmitter regularly after use; never store it wet.
- Keep the Polar Heart Rate Monitor out of extreme cold and heat.
- Do not expose the Polar Heart Rate Monitor to direct sunlight for extended periods.


1. HOW THE POLAR HEART RATE MONITOR WORKS	46
Polar WearLink™ Transmitter	46
How to shorten the wrist band	47
1.1 BUTTONS AND THEIR FUNCTIONS	48
1.2 MODES AND THEIR FUNCTIONS	49
2. SETTINGS	50
2.1 USER INFORMATION SETTINGS	50
Setting User 1 or 2	50
Setting Weight	51
Setting Height	51
Setting Age	51
Setting Sex	51
Setting Activity Level	51
2.2 MEASUREMENT OPTIONS SETTING	52
Turning OwnZone ON/OFF	52
Setting Heart Rate Limits Manually	52
Turning the OwnCal Calorie Counter ON/OFF	53
Turning Polar Fitness Test ON/OFF	54
2.3 SETTING THE WATCH	55
Setting the Alarm	55
Setting the Time of Day	56
Setting the Date	56
3. START THE HEART RATE MEASUREMENT	58
4. OWNZONE	59
4.1 HOW TO DETERMINE YOUR OWNZONE	59
4.2 OWNZONE METHODS	60
4.3 HOW TO USE PREVIOUSLY DETERMINED OWNZONE	61
5. FUNCTIONS DURING HEART RATE MEASUREMENT	62
See Your Heart Rate	62
Illuminate the Display	62
Turn the Sound Signal and the Button Sound ON/OFF	62
Swap the User	62
Start/Stop the Stopwatch	62
Check Exercise Information	63
Check Your Heart Rate Limits	63
See Your Progress Bullets During Exercise	63
Swap the Heart Rate Limits	63
6. STOP THE HEART RATE MEASUREMENT	64
7. RECALLING TRAINING INFORMATION	65
Date	65
Exercise Time (EXE. TIME)	65
Exercise Time in the Target Zone (IN ZONE)	65

Target Heart Rate Limits (LIMITS)	65	10. PRECAUTIONS	73
Average Heart Rate (AVERAGE Hr)	65	10.1 USING A POLAR HEART RATE MONITOR	
Calories Expended During the Exercise (KCAL FAT)	65	IN A WATER ENVIRONMENT	73
Total Calorie Count (Tot. KCAL)	65	10.2 POLAR HEART RATE MONITOR AND INTERFERENCE	73
Total Exercise Time (Tot. TIME)	66	10.3 MINIMIZING POSSIBLE RISKS IN EXERCISING WITH	
Latest OwnIndex	66	HEART RATE MONITOR	74
8. OWNINDEX	67	11. FREQUENTLY ASKED QUESTIONS	75
8.1 FITNESS TEST SETTINGS	67	12. POLAR GLOSSARY	77
Set the Activity Level	67	13. TECHNICAL SPECIFICATIONS	78
Set the Fitness Test Mode On	68	14. LIMITED POLAR INTERNATIONAL GUARANTEE	79
8.2 CARRYING OUT THE TEST	68	15. DISCLAIMER	80
Prepare the Wrist Unit	68	INDEX	81
Relax Yourself for 1–3 Minutes	68		
Start the Test	69		
8.3 FITNESS TEST RESULT	69		
9. CARE AND MAINTENANCE	71		
9.1 TAKING CARE OF YOUR			
POLAR HEART RATE MONITOR	71		
9.2 BATTERIES	71		
9.3 SERVICE	72		

1. HOW THE POLAR HEART RATE MONITOR WORKS


Polar WearLink™ Transmitter:

Electrode areas




Polar WearLink Strap


1. Moisten the electrode areas of the strap under running water and make sure that it is well moistened.



2. Attach the connector to the strap. Position the connector's letter R to the word RIGHT on the strap and snap the fastener.



3. Adjust the strap length to fit snugly and comfortably. Secure the strap around your chest, just below the chest muscles, and snap the second fastener.




4. Check that the wet electrode areas are firmly against your skin and that the Polar logo of the connector is in a central, upright position. It is important that the electrode areas are wet during exercise.

! It is recommended that you wear the transmitter against your bare skin to ensure flawless operation. However, if you wish to wear the transmitter over a shirt, moisten the shirt well under the electrodes.

! See care instructions of the WearLink Transmitter in chapter 9. Care and maintenance, page 71.


How to Shorten the Wrist Band

Wear the wrist unit as a watch. If the wrist band needs to be shortened, visit a local watch or jewellery shop or you can do the shortening yourself by following these instructions.




1. First detach the link cover which is next to the buckle by pressing the spring bar with a sharp-pointed tool.

Buckle
Link cover




2. Take off the link cover.




3. After the link cover is detached, push the spring bar with a pin and pull the spring bar out from the other side.

Spring bar




4. Also detach the second link cover and take off the spring bar.



5. Take off the black extension link. Do not detach the buckle connecting link.

Extension link



6. Attach the buckle connecting link into the extension link and insert the spring bar.

Extension link
Buckle connecting link





7. Put the link cover on top of the extension link so that the spring bar is correctly placed on one side of the link cover. Then press the spring bar from the other side and slide the link cover to its place. You will hear a clicking sound, when the spring bar is placed correctly.

Link cover




8. There are three small holes on both sides of the buckle plate. By shifting the spring bar you can do the final adjustment of the band length.

*Buckle
Adjustment holes*

If you need to shorten the wrist band more repeat the steps 1–7 but remove the link from the opposite side this time so that the buckle will stay in the middle. Please note that you can only detach the extension links with small arrows on their back side.

1.1 BUTTONS AND THEIR FUNCTIONS



Start/Stop/Ok

- Start heart rate measurement
- Start the stopwatch
- Stop the stopwatch
- Activate the displayed mode
- Lock your selection

▲ Up

- Move forward to the next mode
- Select the desired value (up)

▼ Down

- Move backward to the previous mode
- Select the desired value (down)

☰/🔊

- Backlight on
- Signal sound on or off

Reset

- Resets the heart rate monitor, watch settings will return to the default settings.



Briefly pressing a button (approximately 1 second) will do different things than pressing and holding button for a longer period of time (approx. 2 to 5 seconds).

1.2 MODES AND THEIR FUNCTIONS

Your Polar Heart Rate Monitor has four main modes:



1. Time of Day

You can use the Polar Heart Rate Monitor as an ordinary watch with date, weekday indicator and alarm function.



Press and hold START/STOP/OK to return to the Time of day display, wherever you are in the Options setting or File recalling cycle.



Your Polar Heart Rate Monitor starts reminding you at full hours if more than three days has passed since your last recorded exercise. Keep exercising regularly!



2. Exercise

In the Exercise mode you measure the heart rate.

In this mode you can

- see your heart rate reading
- record the exercise data
- use the stopwatch

- determine your OwnZone®
- see your exercising limits
- see the time of day
- see the time you have exercised within your target zone so far
- see the OwnCal® calorie count of your exercise session
- swap heart rate limits
- swap user (before stopwatch is on)



3. File

Polar Heart Rate Monitor saves information on your latest exercise session. In the File mode you can see the following information:

- date of the exercise
- total exercising time
- exercising time spent in your target zone
- heart rate limits used during exercise
- average heart rate
- OwnCal calorie count and fat percentage of the energy expenditure during the exercise session
- total calorie expenditure starting from the moment you have reset the counter
- total exercising time starting from the moment you have reset the counter
- your latest OwnIndex™ based on Polar Fitness Test™



4. Options

In the Options mode you can define the following information:

- User specific information; user 1/2, weight, height, age, sex and activity level for Polar Fitness Test. You have two users and can set the preferences for both user 1 and user 2.
- OwnZone basic, low or high. If you want to determine the values for the target zone manually, set OwnZone OFF.
- OwnCal calorie counter ON/OFF
- Polar Fitness Test ON/OFF
- Watch settings; alarm, 12h/24h mode, time of day, date

2. SETTINGS

In the Options mode you can set

- your personal user information which is a prerequisite for using OwnZone and OwnCal functions
- measurement options which are used during Exercise mode
- watch functions
- Polar Fitness Test function ON or OFF



At any time during the settings you can return to the Time of day display by pressing and holding START/STOP/OK.



The digits run faster if you press and hold the ▲ Up or ▼ Down. Release the button at the desired value.

2.1 USER INFORMATION SETTINGS

Begin with the display showing the time of day.

1. Press ▲ Up or ▼ Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the Options mode.
USER SET is displayed.
3. Press START/STOP/OK to enter the user information loop.

Setting User 1 or 2

If you are in the Time of day display start from **User information**.



You can set two users' preferences. This function is practical e.g. for another family member. Make the settings ready for both users: when you start exercising, just select your personal user id 1 or 2. USER appears, 1 or 2 is flashing.

1. Press ▲ Up or ▼ Down to select USER 1 or USER 2.
2. Press START/STOP/OK to lock your selection.
You have now selected the user id. All settings are made to the selected user.



In File, Measure, Options and Fitness test main modes you can change the user id by pressing and holding ▲ Up.


If you want to disable the user id 2, set the weight of the user 2 to zero.

Continue by setting your weight or return to the Time of day display by pressing and holding START/STOP/OK.

Setting Weight

If you are in the Time of day display start from **User information**. Press START/STOP/OK until LBS or KG appears and the digits of pounds/kilograms are flashing.

1. Press ▲ Up or ▼ Down to adjust your weight.
2. Press START/STOP/OK to lock your selection.
You have now set your weight.


 You can change the weight unit between lbs and kg in watch settings.
12 h time mode -> lbs
24 h time mode -> kg

Continue the user settings or return to the Time of day display by pressing and holding START/STOP/OK.

Setting Height

If you are in the Time of day display start from **User information**. Press START/STOP/OK until FT/INCH or CM appears and the digits of feet/centimeter are flashing.

1. Press ▲ Up or ▼ Down to adjust your height.
2. Press START/STOP/OK to lock your selection.
You have now set your height.

 You can change the height unit between cm and ft/ inch in watch settings.
12 h time mode -> ft
24 h time mode -> cm

Continue by setting your age or return to the Time of day display by pressing and holding START/STOP/OK.

Setting Age

If you are in the Time of day display start from **User information**. Press START/STOP/OK until AGE appears and the digits for years are flashing.

1. Press ▲ Up or ▼ Down to adjust your age.
2. Press START/STOP/OK to lock your selection.
You have now set your age.

Continue by setting your gender or return to the Time of day display by pressing and holding START/STOP/OK.

Setting Sex

If you are in the Time of day display start from **User information**. Press START/STOP/OK until SEX appears and FEMALE/MALE is flashing.

1. Press ▲ Up or ▼ Down to select the right gender.
2. Press START/STOP/OK to lock your selection.
You have now set your gender.

Setting Activity Level

You need to set a suitable activity level for Polar Fitness Test.

If you are in the Time of day display start from **User information**. Press START/STOP/OK until ACTIVITY appears and LOW/MIDDLE/HIGH/TOP is flashing.

1. Press ▲ Up or ▼ Down to select the right activity level.
2. Press START/STOP/OK to lock your selection.

You have now set all necessary user information.
USER SET is displayed.

Choose from the following two alternatives to continue:

1. To set OwnZone function ON or OFF, press ▲ Up.
2. To return to the Time of day display, press ▲ Up or ▼ Down until you come to the Time of day display or simply press and hold START/STOP/OK.

2.2 MEASUREMENT OPTIONS SETTING

Turning OwnZone ON/OFF

Select whether you want to use the OwnZone option. You can select between three different OwnZones. Alternatively, you can set the target zone limits manually yourself.

- OwnZone BASIC for allround exercising
- OwnZone LOW for light intensity exercising
- OwnZone HIGH for moderate intensity exercising
- OFF – The OwnZone function is off and you can set the target zone limits manually.

If you choose to use OwnZone, Polar Heart Rate Monitor asks your user information, if you haven't set them earlier in the USER SET mode.

Begin with the display showing the time of day.

1. Press ▲ Up or ▼ Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the setting loop.
USER SET is displayed.
3. Press ▲ Up once.
OWNZONE OFF or BASIC or LOW or HIGH is displayed.
4. Press START/STOP/OK to enter the setting loop.
OWNZONE is displayed.
OFF or BASIC or LOW or HIGH is flashing.
5. Press ▲ Up or ▼ Down to select a suitable OwnZone for your exercise or OwnZone OFF.
6. Press START/STOP/OK to lock your selection.

Choose from the following alternatives to continue:

1. If you selected OwnZone OFF continue directly from **Setting heart rate limits manually.**
or
2. If you selected OwnZone ON continue from **Turning the OwnCal calorie counter ON/OFF** by pressing ▲ Up once.
or
3. Return to the Time of day display by pressing and holding START/STOP/OK.

Setting Heart Rate Limits Manually

If you choose not to use the OwnZone function, you may manually enter a target zone. When you start the manual target zone settings for the first time:

- Polar Heart Rate Monitor automatically calculates you heart rate limits based on age, if you have set your user information.

- If you haven't done the user settings, default settings 80/160 are shown on the display.

You can start exercising with these limits right away.


Begin with the display showing the time of day.

1. Press ▲ Up or ▼ Down until you enter the Options mode.
2. Press START/STOP/OK to enter the setting loop.
USER SET is displayed.
3. Press ▲ Up once.
OWNZONE OFF is displayed.
4. Press START/STOP/OK to enter the setting loop.
OWNZONE is displayed.
OFF is flashing.
5. Press START/STOP/OK to lock OwnZone OFF.

LIMITS is displayed.

The upper limit is flashing.

6. Press ▲ Up or ▼ Down to adjust the value.

 The digits run faster if you press and hold the ▲ Up or ▼ Down button. Release the button at the desired value.

7. Press START/STOP/OK to lock your selection.
The lower limit reading starts to flash.
8. Press ▲ Up or ▼ Down to adjust the correct value.
9. Press START/STOP/OK to lock your selection.
OWNZONE OFF is displayed.



The highest possible value for the upper limit is 240 and the lowest possible value for the lower limit is 30. The lower limit cannot be set above the upper limit.

Choose from the following two alternatives to continue:

1. To continue setting the OwnCal function ON or OFF, press ▲ Up.
2. To return to the Time of day display, press ▲ Up or ▼ Down until you come to the Time of day display or simply press and hold START/STOP/OK.

Turning the OwnCal Calorie Counter ON/OFF

Begin with the display showing the time of day.

1. Press ▲ Up or ▼ Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the setting loop.
USER SET is displayed.
3. Press ▲ Up or ▼ Down until OWNCAL is displayed.

OWNCAL ON or OFF is displayed.

4. Press START/STOP/OK to enter the setting loop.
OWNCAL is displayed.
ON or OFF is flashing.
5. Press ▲ Up or ▼ Down to select whether you want the calorie expenditure count to be ON or OFF.
6. Press START/STOP/OK to lock your selection.

Choose from the following two alternatives to continue:

1. To continue settings, press ▲ Up.
2. To return to the Time of day display, press ▲ Up or ▼ Down until you come to the Time of day display or simply press and hold START/STOP/OK.

The OwnCal calorie calculation starts when your heart rate reaches 100 bpm. The higher the heart rate, the faster the energy expenditure accumulation. Energy expenditure depends on body weight; the heavier the weight the higher the energy expenditure. Men's calorie expenditure is more intense than women's at the same heart rate.

The energy expenditure measurement with Polar Heart Rate Monitor is most accurate in continuous activities such as walking, cycling and jogging. Polar Heart Rate Monitor shows the energy expenditure in kilocalories (kcal).
1 kilocalorie (kcal or Calorie) = 1000 calories (cal).

When you expend calories, your body actually uses fat and carbohydrates as energy resources. The amount of fat from the total calorie expenditure varies between 10%–60%. The rest (40%–90%) is mostly carbohydrates.

The lower the heart rate, the bigger the fat percentage is of the total energy expenditure. If you want to control weight by exercising, it is recommended to aim at as big a fat percentage of total energy expenditure as possible. Exercise at the lower limit of your OwnZone.

Exercise with OwnZone Low.

Turning Polar Fitness Test ON/OFF

Begin with the display showing the time of day.

1. Press ▲ Up or ▼ Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the setting loop.
USER SET is displayed.
3. Press ▲ Up or ▼ Down until FIT. TEST is displayed.

FIT. TEST ON or OFF is displayed.

4. Press START/STOP/OK to enter the setting loop.
FIT. TEST is displayed.
ON or OFF is flashing.
5. Press ▲ Up or ▼ Down to select whether you want the Polar Fitness Test to be ON or OFF.
6. Press START/STOP/OK to lock your selection.


If you select Polar Fitness Test ON, the Fitness test mode appears in the main mode loop. The Polar Fitness Test is carried out in the Fitness test mode.

Choose from the following two alternatives to continue:

1. To continue settings, press ▲ Up
2. To return to the Time of day display, press ▲ Up or ▼ Down until you come to the Time of day display or simply press and hold START/STOP/OK.

2.3 SETTING THE WATCH





Quick enter to the Watch setting loop: In the Time of day display press and hold the  button for a couple of seconds.

If you want to skip some parts and go directly to set a particular watch function, press START/STOP/OK until the desired watch function is displayed.









To return to the Time of day display at any stage of watch settings, press and hold START/STOP/OK.









The digits run faster if you press and hold the  Up or  Down. Release the button at the desired value.

Setting the Alarm

1. In the Time of day display press  Up or  Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the Options mode. USER SET is displayed.
3. Press  Up or  Down until WATCH SET is displayed.
4. Press START/STOP/OK to start the watch settings. ALARM is displayed. ON or OFF is flashing.
5. Press  Up or  Down to set the alarm ON or OFF.

When the alarm is turned on  appears on the display.

6. Press START/STOP/OK to lock your selection. If you have 12h time mode, AM or PM starts to flash.
7. Press  Up or  Down to select AM or PM.
8. Press START/STOP/OK to lock your selection. The hours start to flash.
9. Press  Up or  Down to adjust the desired hour.
10. Press START/STOP/OK to lock your selection. The minutes start to flash.
11. Press  Up or  Down to adjust the desired minutes.
12. Press START/STOP/OK to lock your selection.

Choose from the following two alternatives to continue:

1. To set the time of day continue in chapter Setting the time of day.
2. To return to the Time of day display, press and hold START/STOP/OK.



When the alarm comes on, it will sound for one minute or you can turn it off by pressing any of the four buttons.

Setting the Time of Day

1. In the Time of day display press ▲ Up or ▼ Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the Options mode.
USER SET is displayed.
3. Press ▲ Up or ▼ Down until WATCH SET is displayed.
4. Press START/STOP/OK to start the watch settings.
ALARM is displayed.
ON or OFF is flashing.
5. Press START/STOP/OK until TIME and the time of day are displayed.

TIME and the time of day are displayed.
12h or 24h is flashing.

6. Press ▲ Up or ▼ Down to select 12h or 24h mode.
7. Press START/STOP/OK to lock your selection.
If you chose 12h time mode, TIME AM or TIME PM is displayed.
8. Press ▲ Up or ▼ Down to select AM or PM.
9. Press START/STOP/OK to lock your selection.
The hours start to flash.
10. Press ▲ Up or ▼ Down to adjust the desired hour.
11. Press START/STOP/OK to lock your selection.
The minutes start to flash.
12. Press ▲ Up or ▼ Down to adjust the desired minutes.
The seconds go to zero.
13. Press START/STOP/OK to lock your selection.
The seconds start running.

Choose from the following two alternatives to continue:

1. To set the date continue in chapter Setting the date.
2. To return to the Time of day display, press and hold START/STOP/OK.

Setting the Date

1. In the Time of day display press ▲ Up or ▼ Down until OPTIONS is displayed.
2. Press START/STOP/OK to enter the Options mode.
USER SET is displayed.
3. Press ▲ Up or ▼ Down until WATCH SET is displayed.
4. Press START/STOP/OK to start the watch settings.
ALARM is displayed.
ON or OFF is flashing.
5. Press START/STOP/OK until DAY or MONTH is displayed.



DAY or MONTH is displayed.

The first digits start to flash.

The date is displayed differently depending on which time mode (12h or 24h) has been chosen.

12h time mode:
MONTH – DAY – YEAR

24h time mode:
DAY – MONTH – YEAR

6. Press ▲ Up or ▼ Down to adjust the desired value for the day/month.
7. Press START/STOP/OK to lock your selection.
MONTH or DAY is displayed.
The middle digits start to flash.
8. Press ▲ Up or ▼ Down to adjust the desired value for the month/day.
9. Press START/STOP/OK to lock your selection.
YEAR is displayed.
The digits for the year start to flash.
10. Press ▲ Up or ▼ Down to adjust the desired value for the year.
11. Press START/STOP/OK to lock your selection.
WATCH SET is displayed.

Choose from the following two alternatives to continue:

1. Press START/STOP/OK to continue from the alarm setting.
2. To return to the Time of day display, press ▲ Up or ▼ Down until you come to the Time of day display or simply press and hold START/STOP/OK.

3. START THE HEART RATE MEASUREMENT

1. Wear the transmitter and the wrist unit as described in the chapter **How the Polar Heart Rate Monitor works**.
2. To make sure that the code search will be successful:
 - Keep the wrist unit within 1 meter/3 feet of your transmitter.
 - Check that you are not near to other people with heart rate monitors, high voltage power lines, televisions, mobile phones, cars, exercise equipment etc.
3. Begin with display showing time of day.
4. Press START/STOP/OK.



♥ MEASURE appears in the display. A heart shaped frame appears. The wrist unit starts searching for heart rate transmission for max. 15 seconds.

♥ Heart rate measurement commences and the heart symbol inside the frame starts flashing. Your heart rate reading (in beats per minute) appears in the display.

If the code search was not successful, the frame around the heart symbol will disappear after 15 seconds. Go back to step 2 to repeat the code search.

If you still do not succeed in code searching, but your heart rate is displayed and the heart symbol without a frame is flashing, you can

start exercising but your wrist unit may become disturbed by other heart rate monitors.

Check that your personal user id is active (USER 1 or USER 2). You can change the user before starting the stopwatch by pressing and holding ▲ Up. To be able to swap the user you need to set preferences for two users.



5. Start the stopwatch by pressing START/STOP/OK.

Stopwatch starts running. Automatic recording begins. If you set OwnZone on: OZ # ----- appears in the display and OwnZone determination begins. If you set heart rate limits manually: EXERCISE appears in the display.

6. You can now start warmup with your Polar Heart Rate Monitor.



Quick start: In the Time of day display, press and hold START/STOP/OK for a couple of seconds.

To stop the heart rate measurement: Press START/STOP/OK to stop the stopwatch. If you want to QUIT the heart rate measurement mode, press ▲ Up or ▼ Down to enter the Time of day display.

4. OWNZONE

Polar Heart Rate Monitor is able to determine the individual training heart rate zone for you. This exercise zone is called your OwnZone (OZ).

For most adults, OwnZone corresponds to 65–85 % of the maximum heart rate. OwnZone is suitable for cardiovascular training either for health-related or fitness improvement purposes.

OwnZone works well for those whose goal is weight control via physical activity and exercise, because the OwnZone corresponds to intensity from "light to moderate" to "moderate". Exercise of this intensity level can be conducted even for extended periods safely and often yields the best results for weight loss.

The OwnZone can be found during a warmup period of an exercise session lasting 1–5 minutes.

4.1 HOW TO DETERMINE YOUR OWNZONE

OwnZone is easily determined by walking and jogging as instructed below. However, you are able to find your OwnZone with other exercise modes as well, e.g. exercising with a stationary bike. The principle is that you start exercising very slowly with light intensity and gradually increase your heart rate.

1. Check that you have set your personal user information correctly.

2. Because OwnZone determination needs every heart beat to be picked up properly, make sure that the transmitter belt is tight enough and the electrodes well moistened.
3. Check that you have set the OwnZone function on.
4. Start the heart rate measurement as described earlier.
5. Start the stopwatch.
OZ # _ _ _ _ appears in the display.
Stopwatch starts running.

6. There are five stages in finding your OwnZone. OwnZone appears typically on the 3rd or 4th stage (in 3–4 minutes). Keep your heart rate below 100 bpm during the first stage. After each stage you should increase your heart rate 10–20 bpm. The display shows which stage you are in:

- | | |
|-------------|---------------------------------------|
| # _ _ _ _ _ | Walking at very slow pace for 1 min. |
| # # _ _ _ _ | Walking at brisk pace for 1 min. |
| # # # _ _ _ | Walking at very brisk pace for 1 min. |
| # # # # _ _ | Jogging at slow pace for 1 min. |
| # # # # # | Running for 1 min. |

7. At some point during the stages 1–5 you hear two beeps. This indicates that Polar Heart Rate Monitor has determined your exercising limits and you can continue your exercising.



The determined heart rate limits are displayed for a few seconds. EXERCISE appears in the display.

! You can recall the limits onto the display by bringing the wrist unit close to the transmitter.

You can turn the OwnZone timing signals off by pressing and holding \Rightarrow (right arrow) for two seconds before starting the OwnZone determination.

Reset the stopwatch:

If needed you can reset the stopwatch to zero and restart OwnZone determination:

1. Stop the OwnZone determination by pressing START/STOP/OK. Heart rate limits for your previously determined OwnZone appear on the display.
2. Stop the stopwatch by pressing START/STOP/OK. QUIT? appears in the display.
3. To reset the stopwatch press and hold \blacktriangledown Down. The stopwatch will start blinking and RESET? text appears in the display.
To reset the stopwatch keep pressing \blacktriangledown Down.
4. To restart the stopwatch press START/STOP/OK.

4.2 OWNZONE METHODS

Polar Heart Rate Monitor uses multiple methods for determining your OwnZone. The primary method is based on measuring and analysing the heart rate variability during graded exercise (OZ.V).

If your heart rate exceeds the safety limit during the determination, Polar Heart Rate Monitor picks up your latest OwnZone definition based on heart rate variability (OZ.L). Alternatively, if there is no previous determination in memory, OwnZone limits, which are based on your age predicted maximum heart rate (OZ.A), will be used. These methods ensure that the safest possible exercise zone is determined for you.

The method used in defining the OwnZone is indicated after the two beeps as follows:

1. OZ.V (OwnZone Variability) and the currently determined OwnZone limits appear on the illuminated display.
2. OZ.L (OwnZone Latest) and your previously determined OwnZone limits appear on the illuminated display.
3. OZ.A (OwnZone Age) and your age-predicted heart rate limits appear on the illuminated display.

You will also see which OwnZone level you have selected to use:
BASIC (OwnZone Basic for allround exercising)
HIGH (OwnZone High for exercising at moderate intensity)
LOW (OwnZone Low for exercising at light to moderate intensity)

4.3 HOW TO USE PREVIOUSLY DETERMINED OWNZONE

Every time you start the stopwatch and have the OwnZone option on, the Polar Heart Rate Monitor starts finding your OwnZone. You can determine your OwnZone for every exercise session separately. At least it should always be performed when changing exercise environment, exercise mode or if you have abnormal feelings about your physiological state when starting the exercise.

However, you can choose to exercise with the latest OwnZone limits in case you are not willing to start the process redefining your OwnZone.

1. Start the heart rate measurement and stopwatch as usual.
OZ # ---- appears on the display.
2. Press START/STOP/OK. Polar Heart Rate Monitor starts using your latest determined OwnZone.
OZ.L and the values for your previously determined OwnZone appear momentarily on the illuminated display.

In case there is no previously determined OwnZone in memory, Polar Heart Rate Monitor starts using your age predicted exercising limits.

OZ.A and your age-predicted heart rate limits appear on the illuminated display.


5. FUNCTIONS DURING HEART RATE MEASUREMENT

MEASURE



See Your Heart Rate

You can see your heart rate in beats per minute. A flashing heart symbol indicates an ongoing heart rate measurement. The inner heart symbol flashes at the pace of your heart. The outline indicates that the heart rate reception is coded.


 If there is - - reading on the display, there has not been heart rate reception for at least 5 minutes. Bring the wrist unit up to your chest near the transmitter's Polar logo. The wrist unit starts looking for the heart rate again.

Illuminate the Display


Press 

Backlight is activated for a few seconds.

Turn the Sound Signal and the Button Sound ON/OFF

When  appears in the display, you will hear an alarm sound with each heart beat if you go above the upper limit or below the lower limit of your Target Zone. You will hear a beep when pressing the buttons.

The sound signal will be on in OwnZone determination. The sound signal will also alarm you at the end of the Fitness test.

Press and hold  for 2 seconds to turn the sound signal and the button sound on and off.

-)) Sound signal is on.



A flashing heart rate reading also informs you of being outside your target zone.

The sound signal setting is saved i.e. the next time you start heart rate measurement the previous setting will be used.

Swap the User

You can swap the user in the beginning of an exercise when the stopwatch is not running yet.

To be able to swap the user you need to set preferences for two users.

Press and hold 

 or  is displayed to indicate which one of the user settings are activated.

The following exercise will be recorded in the selected user's file.

Start/Stop the Stopwatch

Press START/STOP/OK to start/stop the stopwatch.

Exercise information will be stored in an exercise file only when the stopwatch is running.

You can pause the stopwatch by pressing START/STOP/OK. Repeating it will restart the stopwatch and the above calculations.

To reset the stopwatch to zero:

1. Stop the stopwatch.
2. Press and hold ▼ Down.

Check Exercise Information

When the stopwatch is running press ▲ Up or ▼ Down to see the following information during the exercise:



Time of day (TIME)



Time spent in target zone (IN ZONE)



Calorie expenditure (OWNCAL)

The calorie counting feature must be set on to be able to see the OWNCAL display.



Stopwatch (EXERCISE)

Check Your Heart Rate Limits

While exercising you don't have to press any button to check the heart rate limits you are using. Just bring the wrist unit up to your chest near the transmitter's Polar logo.

The current exercising zone values and backlight will be on for a few seconds.

See Your Progress Bullets During Exercise



A bullet appears for every 10 minutes exercise in Target Zone.

Swap the Heart Rate Limits

You can swap heart rate limits during the exercise while the stopwatch is running and you have set OwnZone on. You can select between OwnZone Basic, OwnZone High, OwnZone Low and manually set limits.

Press and hold ▲ Up or ▼ Down.



If you swap heart rate limits during exercise, the latest ones in use will be saved as measurement option settings and your next exercise will start with these limits.

6. STOP THE HEART RATE MEASUREMENT

1. Press START/STOP/OK to stop the stopwatch.
The stopwatch stops.

QUIT?  appears in the display.

Heart rate measurement without exercising limits or calorie calculation continues.

2. Press ▲ Up or ▼ Down to leave the Measure mode.
Heart rate measurement stops.



If you press and hold ▼ Down for 1.5 seconds, the stopwatch will start blinking and RESET? text appears in the display. You can restart the stopwatch by pressing START/STOP/OK. Or, if you want to reset the stopwatch keep pressing ▼ Down. The stopwatch will be reset and the current exercise information will be deleted.



The Heart Rate Monitor wrist unit turns itself into Time of day mode within 5 minutes in case you forget to quit the heart rate measuring mode after having stopped the stopwatch and taken the transmitter off your chest.

7. RECALLING TRAINING INFORMATION

Your latest training information is in the memory of the Polar Heart Rate Monitor until the next time you start the stopwatch. Then your previous training information will be replaced by the new one. Polar Heart Rate Monitor saves only one training file at a time. The training information is saved, if the stopwatch has been on more than one minute.

1. In the Time of day display, press ▲ Up or ▼ Down until you enter the FILE mode.
FILE is displayed.
2. Press START/STOP/OK to enter the training information recalling loop.
3. Press ▲ Up or ▼ Down to recall the training information.

Date

The date of your last exercise session.

Exercise Time (EXE. TIME)

The exercise time is the time which you have exercised with the stopwatch running.

Exercise Time in the Target Zone (IN ZONE)

The exercise time, which you have spent within your OwnZone or your self-defined target zone.

Target Heart Rate Limits (LIMITS)

You will see the target heart rate limits of your exercise.

Average Heart Rate (AVERAGE Hr)

You will see the average heart rate of your exercise session.



Calories Expended During the Exercise (KCAL/ FAT)

You will see how many kilocalories you have expended during the exercise and the fat percentage of the total energy expenditure.

Total Calorie Count (Tot. KCAL)

The total calorie counter is a function, which memorises information for more than one exercise session. You will have a count for cumulative calories expended during several exercise sessions. This function allows you to follow the calorie count as a gauge of effective exercise e.g. during one week.

You can subtract the total calorie amount.

The kilocalories are subtracted by 10-kcal decrement.

1. Start with the display showing Tot. KCAL.
2. Press START/STOP/OK to enter the Tot. KCAL mode.
The amount of total kilocalories starts flashing.
3. You can now subtract the calories by pressing ▼ Down.
4. Press START/STOP/OK to lock the desired amount.



Once you have subtracted the Tot. KCAL value, you cannot return it.



This function is useful in weight management, if you want to follow the combined effects of exercise and nutrition.

To reset the total calorie expenditure count to zero:

1. Start with the display showing Tot. KCAL.
2. Press START/STOP/OK to enter the Tot. KCAL mode.
The amount of total calories starts flashing.
3. Press and hold ▼ Down for a couple of seconds.
RESET starts to flash.
0 appears.

If you release the ▼ Down button, when RESET is flashing, the value is not reset.

Total Exercise Time (Tot. TIME)

The cumulative exercise time counter is a function, which memorizes information for more than one exercise session. You will have a count for the cumulative time of several exercise sessions. This function can be used as a gauge of exercise time e.g. during one week.

To reset the cumulative exercise time count to zero:

1. Start with the display showing Tot. TIME.
2. Press START/STOP/OK to enter the Tot. TIME mode.
The amount of cumulative exercise time starts flashing.
3. Press and hold ▼ Down for a couple of seconds.
RESET starts to flash.
00:0000 appears.

If you release the ▼ Down button, when RESET is flashing, the value is not reset.

Latest OwnIndex

You can see your OwnIndex based on your latest Polar Fitness Test and the date when it was carried out.

8. OWNINDEX

OwnIndex is a score which is comparable to VO_{2max} , a commonly used descriptor of aerobic fitness. OwnIndex is a result of the Polar Fitness Test, which is an easy, safe and quick way to estimate maximal aerobic power. The OwnIndex is based on heart rate, heart rate variability at rest, body structure and self-assessed physical activity.

Cardiovascular or aerobic fitness relates to how well your cardiovascular system works to transport oxygen to your body. The better the cardiovascular fitness, the stronger and more efficient your heart is.

Having good cardiovascular fitness has many health benefits. For example, it decreases your risk of cardiovascular diseases, stroke and high blood pressure.

If you want to improve your cardiovascular fitness, it takes a minimum of 6 weeks of regular training to see a noticeable change. Less fit individuals see progress even more rapidly and for more active individuals more time is needed.

Cardiovascular fitness is best improved by exercise types which use large muscle groups. Such activities include e.g. walking, jogging, swimming, rowing, skating, cross-country skiing and cycling.

To be able to follow your fitness progress start with measuring your OwnIndex a couple of times during the first two weeks to get a reliable baseline value. Thereafter repeat the test approximately once a month.

Polar Fitness Test is targeted to healthy adults.

8.1 FITNESS TEST SETTINGS

To be able to carry out the Polar Fitness Test you need to do the following settings:

- Select your long-term physical activity level.
- Set the Fitness test mode on.

Set the Activity Level

Begin with the display showing the time of day.

1. Press ▲ Up or ▼ Down until you enter the Options mode.
2. Press START/STOP/OK to enter the setting loop.
USER SET is displayed.
3. Press START/STOP/OK to enter the user information loop.
4. Press START/STOP/OK until ACTIVITY appears.
ACTIVITY is displayed and LOW/MIDDLE/HIGH/TOP is flashing.
5. Press ▲ Up or ▼ Down to select the right activity level.
6. Press START/STOP/OK to lock your selection.
USER SET is displayed.
7. To return to the Time of day display, press and hold START/STOP/OK.

Activity Level

Assess your long-term physical activity level. Do not change your activity level description if your regular exercise habits have changed less than six months ago.

- Low** Exercise is not a regular part of your life. You walk only for pleasure or occasionally exercise sufficiently to cause heavy breathing or perspiration.

- Middle** You participate regularly in recreational sports. For example, you run 5 miles/8km a week or spend 30–60 minutes a week in comparable physical activity. Or, your work requires modest physical activity.
- High** Exercise is an inseparable part of your life. You exercise regularly at least 3 times a week with heavy paced intensity. Exercise causes clear breathlessness and sweating. For example, you run regularly more than 10 miles/16 km a week or spend more than 1.5–2 hours in comparable physical activity.
- Top** You participate regularly in heavy physical exercise at least 5 times a week. E.g. you exercise to improve performance for competitive purposes.

Set the Fitness Test Mode On

1. In the Time of day display press ▲ Up or ▼ Down until you enter the Options mode.
OPTIONS is displayed.
2. Press START/STOP/OK to enter the setting loop.
USER SET is displayed.
3. Press ▲ Up or ▼ Down until you come to the Fitness test setting mode.
FIT. TEST ON/OFF is displayed.
4. Press START/STOP/OK to enter the Fitness test setting loop.
FIT. TEST is displayed.
ON or OFF is flashing.
5. Press ▲ Up or ▼ Down to select the test on.

6. Press START/STOP/OK to lock your selection.
7. To return to the Time of day display, press and hold START/STOP/OK.

8.2 CARRYING OUT THE TEST

To get reliable test results, there are the following basic requirements:

- You should be relaxed and calm.
- The test can take place anywhere – home, office, health club – as long as the testing environment is peaceful. No disturbing noises (e.g. television, radio or telephone), no other people talking to you.
- Keep the testing place, time of day and environment the same every time the test is repeated.
- Avoid eating a heavy meal or smoking 2–3 hours prior to testing.
- Avoid heavy physical effort, alcoholic beverages or pharmacological stimulants on the test day and the day before.

Prepare the Wrist Unit

1. In the Time of day display press ▲ Up or ▼ Down until you come to the Fitness test mode.
FIT. TEST is displayed.



Relax Yourself for 1–3 Minutes

2. Lay down, relax your muscles and mind.

Start the Test

- Press START/STOP/OK to start the test.
The wrist unit starts searching for your heart rate.
--- TEST ON and your heart rate are displayed.
The test begins.



If you want to have the sound signal on during the Fitness test, check there is  in the display. Press and hold  to have sound signal on. The sound signal will alarm you at the end of the test.

- It takes 3–5 minutes to carry out the test. Lie relaxed, don't raise arms or legs, avoid body movements during the test. Don't communicate with other people.
Your current OwnIndex interpretation and the testing date are displayed.
- To return to the Time of day display after the test, press and hold START/STOP/OK.

You can stop the test anytime during the test by pressing START/STOP/OK. The latest OwnIndex will not be replaced.


FAILED TEST is displayed for a moment.

Your latest OwnIndex interpretation and the test date are displayed.

If the wrist unit does not succeed to receive your heart rate at the beginning or during the test, the test will fail. Check that the transmitter electrodes are wet and the elastic strap is snug enough. Restart the procedure.

8.3 FITNESS TEST RESULT

The test result, your OwnIndex, is most meaningful when comparing your individual values and changes in it. OwnIndex can also be interpreted according to your gender and age. Locate your OwnIndex on the following table to find out your current cardiovascular fitness classification, compared to those who are of same age and gender.



- 1–3 Increase exercise for health benefits and fitness.
- 4 Maintain current exercise habits to ensure better health. Exercise increase recommended for fitness improvements.
- 5–7 Maintain current exercise habits for good health and fitness.

FITNESS CLASSES	Age	1	2	3	4	5	6	7
	(years)	Very poor	Poor	Fair	Average	Good	Very good	Excellent
MEN	20-24	< 32	32-37	38-43	44-50	51-56	57-62	> 62
	25-29	< 31	31-35	36-42	43-48	49-53	54-59	> 59
	30-34	< 29	29-34	35-40	41-45	46-51	52-56	> 56
	35-39	< 28	28-32	33-38	39-43	44-48	49-54	> 54
	40-44	< 26	26-31	32-35	36-41	42-46	47-51	> 51
	45-49	< 25	25-29	30-34	35-39	40-43	44-48	> 48
	50-54	< 24	24-27	28-32	33-36	37-41	42-46	> 46
	55-59	< 22	22-26	27-30	31-34	35-39	40-43	> 43
60-65	< 21	21-24	25-28	29-32	33-36	37-40	> 40	
WOMEN	20-24	< 27	27-31	32-36	37-41	42-46	47-51	> 51
	25-29	< 26	26-30	31-35	36-40	41-44	45-49	> 49
	30-34	< 25	25-29	30-33	34-37	38-42	43-46	> 46
	35-39	< 24	24-27	28-31	32-35	36-40	41-44	> 44
	40-44	< 22	22-25	26-29	30-33	34-37	38-41	> 41
	45-49	< 21	21-23	24-27	28-31	32-35	36-38	> 38
	50-54	< 19	19-22	23-25	26-29	30-32	33-36	> 36
	55-59	< 18	18-20	21-23	24-27	28-30	31-33	> 33
60-65	< 16	16-18	19-21	22-24	25-27	28-30	> 30	

The classification is based on literature review of 62 studies where VO_{2max} was measured directly in healthy adult subjects in the USA, Canada and 7 European countries. Reference: Shvartz E, Reibold RC: Aerobic fitness norms for males and females aged 6 to 75 years: a review. *Aviat Space Environ Med*; 61:3-11, 1990.

9. CARE AND MAINTENANCE

Your Polar Heart Rate Monitor is a high-tech instrument of superior design and workmanship and should be treated with care. The suggestions below will help you fulfill the warranty obligations and enjoy this product for many years to come.

9.1 TAKING CARE OF YOUR POLAR HEART RATE MONITOR

Connector

- Clean the connector regularly after use with a mild soap and water solution.
- Never use alcohol or any abrasive material such as steel wool or cleaning chemicals on any part of the transmitter.
- Do not dry the connector in any other way than with a towel. Rough handling may damage it.
- Never put the connector in a washing machine or a drier!

Strap

- You can wash the strap in a washing machine at 104 ° F/ +40 °C. We recommend that you use a washing pouch.
- We recommend that you wash the strap after each time you have used it in pool water with a high chlorine content.
- Do not spin-dry the strap!
- Do not put the strap in a dryer!
- Do not iron the strap!

Maintenance

- Store your transmitter in a cool and dry place. Do not store it in any kind of non-breathing material, such as a plastic bag or a

sports bag, if it is wet. Sweat and moisture can keep its electrodes wet and the transmitter activated, which shortens the battery life.

- Keep your transmitter out of extreme cold (below 14 °F/ -10 °C) and heat (above 122 °F/ +50 °C).
- Do not expose transmitter to direct sunlight for extended periods, such as by leaving it in a car.

9.2 BATTERIES


Transmitter

The estimated average battery life of Polar WearLink Transmitter is 2 years (1h/day, 7 days/week). If your transmitter stops working, it might be due to the battery running out. We recommend that you have the battery replaced by an authorized Polar Service Center only.

You can, however, change the battery of the Polar WearLink Transmitter yourself by following the instructions below. To change the battery, you need a coin, sealing ring (some included in the product package) and battery (CR2025).

1. Open the connector's battery cover by turning it counter-clockwise with the coin from the CLOSE position to the OPEN position. The cover pops out a little from its niche.
2. Remove the battery cover from the connector, lift the battery from the niche with your fingers and replace it with a new battery.
3. Remove the sealing ring of the battery cover and replace it with a new ring.

4. Place the negative (-) side of the battery against the connector and the positive (+) side against the cover.
5. Put the battery in the connector's niche and put the cover on so that the arrow points to the OPEN position. Make sure that the cover's sealing ring is placed correctly in its groove.
6. Gently press the cover deep enough so that its exterior surface is on the same level as the connector's surface.
7. Turn the cover with the coin clockwise so that the arrow on the cover turns from OPEN to CLOSE. Make sure that the cover is closed properly!

 In order to ensure the maximum life span of the connector's battery cover, open it only when you need to change the battery. You can get the sealing ring/battery kits from well-equipped Polar retailers and authorized Polar Services. In the USA and Canada: Additional sealing rings are available from an authorized Polar Service Center.


Polar Wrist Unit

The estimated average battery life of the Polar wrist unit is 20 months in normal use (1 h / day, 7 days a week). Note also that excessive use of the backlight and the alarm signals will use the battery more rapidly.

Do not open the Polar wrist unit yourself. To ensure the water resistance properties and the use of qualified components, the wrist unit battery should be replaced only by an authorized Polar Service Center. At the same time a full periodic check of the Polar Heart Rate Monitor will be done.

9.3 SERVICE

We recommend that you have all service done by an authorized Polar Service Center only.

 Water resistance cannot be guaranteed after unauthorized service.

10. PRECAUTIONS

10.1 USING A POLAR HEART RATE MONITOR IN A WATER ENVIRONMENT

Your Polar Heart Rate Monitor is water resistant to 330 feet/100 meters. To maintain the water resistance, do not operate the buttons of the wrist unit under water.

Heart rate measurement in a water environment is technically demanding for the following reasons:

- Pool water with a high chlorine content and seawater are very conductive. The electrodes of a transmitter may short-circuit, which prevents ECG signals from being detected by the transmitter.
- Jumping into the water or strenuous muscle movement during competitive swimming may cause water resistance that shifts the transmitter on the body to a location where it is not possible to pick up ECG signal.
- The ECG signal strength varies depending on the individual's tissue composition and the percentage of people who have problems in heart rate measuring is considerably higher in a water environment than in other use.

10.2 POLAR HEART RATE MONITOR AND INTERFERENCE

Electromagnetic Interference

Disturbances may occur near high voltage power lines, traffic lights, overhead lines of electric railways, electric bus lines or trams, televisions, car motors, bike computers, some motor driven

exercise equipment, cellular phones or when you walk through electric security gates.

Crosstalk

The Polar wrist unit in non-coded mode of operation picks up transmitter signals within 3 feet/1 meter. Non-coded signals from more than one transmitter picked up simultaneously can cause an incorrect readout.

Exercise Equipment

Several pieces of exercise equipment with electronic or electrical components such as LED displays, motors and electrical brakes may cause interfering stray signals. To try to tackle these problems, relocate the Polar wrist unit as follows:

1. Remove the transmitter from your chest and use the exercise equipment as you would normally.
2. Move the wrist unit around until you find an area in which it displays no stray reading or flashing of the heart symbol. Interference is often worst right in front of the display panel of the equipment, while the left or right side of the display is relatively free of disturbance.
3. Put the transmitter back on the chest and keep the wrist unit in this interference-free area as far as it is possible.
4. If the Polar Heart Rate Monitor still does not work with the exercise equipment, this piece of equipment may be electrically too noisy for wireless heart rate measurement.

10.3 MINIMIZING POSSIBLE RISKS IN EXERCISING WITH A HEART RATE MONITOR

Exercise may include some risk, especially for those who have been sedentary.

Before starting a regular exercise program you are recommended to answer to the following questions for health status checking. If the answer is yes to any of the questions, we recommend to consult a doctor before starting an exercise program.

- Have you not exercised for the past 5 years?
- Do you have high blood pressure?
- Do you have high blood cholesterol?
- Do you have symptoms of any disease?
- Are you taking any blood pressure or heart medication?
- Do you have a history of breathing problems?
- Are you recovering from a serious illness or medical treatment?
- Do you use a pacemaker or another implanted electronic device?
- Do you smoke?
- Are you pregnant?

Note that in addition to exercise intensity, medications for heart, blood pressure, psyche, asthma, breathing etc. as well as some energy drinks, alcohol and nicotine, can affect heart rate.

It is important to be sensitive to your body's reactions during exercise. If you feel unexpected pain or excessive fatigue when exercising, it is recommended to stop the exercise or continue at a lighter intensity.

Notice to persons with pacemakers, defibrillator or other implanted electronic device. Persons who have a pacemaker use the Polar Heart Rate Monitor at their own risk. Before starting use, we always recommend an exercise test under doctor's supervision. The test is to ensure the safety and reliability of the simultaneous use of the pacemaker and the heart rate monitor.

If you are allergic to any substance that comes into contact with the skin or if you suspect an allergic reaction due to using the product, check the listed materials in the "Technical Specifications" chapter. To avoid any skin reaction risk with the transmitter, wear it over a shirt. However, moisten the shirt well under the electrodes to ensure flawless operation.



The combined impact of moisture and intense abrasion may cause a black color, which might stain light-colored clothes especially, to come off the transmitter's surface.



If you use insect repellent on your skin, you must ensure that it does not come into contact with the transmitter!

11. FREQUENTLY ASKED QUESTIONS

What should I do if...

...OwnZone determination only seldom gives me heart rate limits which are based on heart rate variability (OZ.V)?

1. Check that the elastic strap is snug enough.
2. Make sure that the electrodes of the transmitter are moistened against your skin.
3. Put on the transmitter 15–20 minutes prior to the start of your exercise to ensure good adaptation between the electrodes and your skin.
4. Start the first stage of determining your OwnZone sufficiently slowly. If you obtain the heart rate limits based on age (OZ.A) within the first 2 minutes you have raised your heart rate too quickly. The heart rate reading should not rise over 20 beats/minute from the starting level. During the first stage you can, for example, put on your sports clothing.
5. If you reach the OwnZone limits based on age (OZ.A) or the previously determined OwnZone limits based on your heart rate variability (OZ.L) only after 5 minutes from starting specification you have not raised your heart rate sufficiently rapidly. Raise your heart rate more quickly during the OwnZone determination.
6. If you have always used walking and running to determine your OwnZone, walking uphill, cycling or using an ergometer may work better for you.
7. Stress or heavy mental loading may reduce your heart rate variability so that the needed point where the heart rate variability vanishes cannot be measured.

Relax for 5–10 minutes before exercising by sitting down and decreasing your heart rate.

8. For a small percentage of people the heart rate variability will not totally disappear with increasing workload. In this case the age-predicted limits are the best estimation for the exercising zone.

...the OwnZone limits based on heart rate variability (OZ.V) vary considerably?

Heart rate variability reacts sensitively to all changes taking place in the body and environment. The OZ.V heart rate limits are affected by, amongst others, your daily state of fitness, the time of OwnZone determination, the type of exercise and state of mind. The lower limit of the OZ.V heart rate limits can change for the same person by as much as 30 beats.

However it is important that the feeling of exercise carried out in the OwnZone is light or moderate. For example if you raise your heart rate too rapidly after the first stage, you may get too high OwnZone heart rate limits. If the OwnZone limits seem to be too high we recommend that the OwnZone determination is done again.

...there is no heart rate reading (- -)?

1. Check that the electrodes of the transmitter are moistened and that you are wearing it as instructed.
2. Check that you have kept the transmitter clean.
3. Check that there are no sources of electromagnetic radiation in the close vicinity of the Polar wrist unit, such as TV sets, mobile phones, CRT monitors etc.

4. Have you had a cardiac event which may have altered your ECG waveform? In this case consult your physician.

...heart symbol flashes irregularly?

1. Check that your Polar wrist unit is inside the transmission range and not further than 3 feet/1 meter from the Polar Transmitter you are wearing.
2. Check that the elastic strap has not become loose during exercise.
3. Make sure that the electrodes of the transmitter are moistened.
4. Make sure that there is no other heart rate transmitter within the reception range (3 feet/1 meter).
5. Cardiac arrhythmia may cause irregular readings. In this case consult your physician.

...other person with heart rate monitor is causing interference?

Your training partner may have coincidentally activated exactly the same code as you. In this case keep distance to your training partner or proceed as follows:

1. Stop the stopwatch and heart rate measurement.
2. Take the transmitter off your chest for 10 seconds.
3. Put the transmitter back and start the heart rate measurement again.

...heart rate reading becomes erratic or extremely high?

You may have come within range of strong electromagnetic signals which cause erratic readings on the wrist unit display. Check your surroundings and move further away from the source of disturbance.

...there are no reactions to any buttons?

Reset the Polar Heart Rate Monitor. Resetting clears watch settings which will return to the default settings. User information and measurement options settings will be saved.

1. Press the RESET button with a pen tip.
Display becomes full of digits.
2. Press any of the four buttons once.
The Polar Heart Rate Monitor enters the Time of day display.

If you don't press any button after reset within one minute the wrist unit goes to the battery saving mode.

...the display is blank or fading?

The Polar Heart Rate Monitor is in battery saving mode when it leaves the factory. Wake up by pressing START/STOP/OK button twice. The Polar Heart Rate Monitor enters the Time of day display.

Usually the first sign of a discharged battery is fading of the digits when using the backlight. Have the batteries checked.

...the battery of the wrist unit must be replaced?

We recommend having all service done by an authorized Polar Service Center. This is especially necessary to ensure the warranty remains valid and is not affected by incorrect repair procedures done by an unauthorized agent. Polar Service will test your wrist unit for water resistance after battery replacement and make a full periodic check of the complete Polar Heart Rate Monitor.

12. POLAR GLOSSARY



Indicates coded transmission of heart rate.

The inner heart flashes to the beat of your heart.

•)) In the Time of day mode indicates that the alarm is on. In the measure mode indicates that the sound signal will alarm if you are exercising outside your Target Zone.

#----- Progress bullets in the measure mode.

One bullet = 10 minutes in target zone.

AGE: Insert your age in User settings.

ALARM: Indicates the alarm clock in watch settings.

AM/PM: 12 hour mode of time.

AVERAGE Hr: Average heart rate of the exercise.

Coded heart rate transmission: Polar coded transmitter automatically locks in a code to transmit your heart rate to your wrist unit. In coded transmission the wrist unit accepts heart rate from your Polar coded transmitter only. Coding significantly reduces crosstalk caused by other users of heart rate monitors. However, it does not necessarily reduce all environmental interference.

Electrodes: The electrode areas are on the backside of the transmitter.

EXERCISE: Heart rate measurement mode, stopwatch is on and the exercise information is recorded.

EXE. TIME: Total exercise time.

FAT: Indicates the fat percentage of the total calorie expenditure.

FILE: File mode to recall training information.

Heart rate: The number of heart beats per minute (bpm)

Heart rate variability: Fluctuations of interbeat intervals. Heart rate variability of a healthy subject is in general large in resting conditions and during light exercise. Heart rate variability gradually disappears when exercise becomes more intense and heart rate increases.

HR_{max}: Maximum heart rate.

IN ZONE: Exercising time in Target Zone.

KCAL: Indicates the energy expenditure of your exercise session.

LIMITS: Manually set target zone limits are determined as percentages of age predicted or measured maximum heart rate. The percentages are converted into heart rate beats per minute.

MEASURE: Heart rate measurement mode, stopwatch is off.

OPTIONS: Setting mode for user information, heart rate measurement options and watch functions.

OWNCAL: Calorie expenditure during exercise.

OWNINDEX: OwnIndex is a score which is comparable to VO_{2max} . OwnIndex is a result of the Polar Fitness Test.

OWNZONE, OZ.: Target zone which Polar Heart Rate Monitor determines for you. OwnZone is an individual guide for light to moderate exercise intensity.

OZ # ----- Indicates OwnZone determination procedure.

OZ.A: Indicates OwnZone limits which are determined according to your age.

OZ.L: Your latest OwnZone limits.

OZ.V: Indicates currently measured OwnZone limits which are based on your heart rate variability.

SEX: Indicates the gender (Male or Female) in User settings.

Target zone, TZ: The area between the upper and lower target zone limits. Can be defined automatically by Polar Heart Rate Monitor or manually by yourself.

TIME: Indicates the time of day in the watch settings.

Time of day mode: Display showing time of day.

Tot. KCAL: Cumulative total calorie expenditure of your exercise.

Tot. TIME: Cumulative total exercise time of your exercise sessions.

USER SET: Indicates the setting loop for your personal information.

13. TECHNICAL SPECIFICATIONS

Polar Heart Rate Monitor is designed to indicate the level of physiological strain and intensity in sports and exercise. No other use is intended or implied. Heart rate is displayed as number of heart beats per minute (bpm).

Polar Connector

Battery type:	CR 2025
Battery sealing ring:	O-ring 20.0 x 1.0 Material FPM
Battery life:	Average 1300 hours of use
Operating temperature:	14 °F to 122 °F / -10 °C to +50 °C
Material:	Polyamide
Water resistance:	To 100 feet/30 meter

Strap

Material:	Polyurethane, polyamide, nylon, polyester and elasthane
-----------	---

Wrist Unit

Battery type:	CR 2032
Battery life:	Average 20 months (1 h / day, 7 days / week)
Operating temperature:	14 °F to 122 °F/-10 °C to +50 °C
Water resistance:	to 330 feet/100 meters
Accuracy of formula for age based OwnZone heart rate limits (OZ.A):	± 1bpm
Accuracy of heart rate measurement:	± 1% or ± 1 bpm, whichever larger, definition applies to steady state conditions

Back Case, Extend Link Cover's and wrist strap buckle:

Stainless steel complying with the EU Directive 94/27/EU and its amendment 1999/C 205/05 on the release of nickel from products intended to come into direct and prolonged contact with the skin.

Default Settings

Time of day	24h
Alarm	OFF
User ID	1
Weight	00
Height	00
Age	0
Sex	Male
Activity	Low
OwnZone	OFF
Upper Limit (manual)	160
Lower Limit (manual)	80
OwnCal	OFF
Polar Fitness Test	OFF

Limit Values

Chronometer	23 h 59 min 59 s
Heart rate limits	30–240 bpm
Kcal	99999 kcal
Tot. kcal	999999 kcal
Tot. time	9999 h 59 min 59 s

14. LIMITED POLAR INTERNATIONAL GUARANTEE

- This limited Polar international guarantee is issued by Polar Electro Inc. for the consumers who have purchased this product in the USA or Canada. This limited Polar international guarantee is issued by Polar Electro Oy for the consumers who have purchased this product in other countries.
- Polar Electro Inc. / Polar Electro Oy guarantees to the original consumer/ purchaser of this product that the product will be free from defects in material or workmanship for two years from the date of purchase.
- **Please keep the receipt or International Guarantee Card, which is your proof of purchase!**
- The guarantee does not cover the battery, damage due to misuse, abuse, accidents or non-compliance with the precautions; improper maintenance, commercial use, cracked or broken cases.
- Guarantee does not cover any damage/s, losses, costs or expenses, direct, indirect or incidental, consequential or special, arising out of, or related to the product. During the guarantee period the product will be either repaired or replaced at an authorized service center free of charge.
- This guarantee does not affect the consumer's statutory rights under applicable national or state laws in force, or the consumer's rights against the dealer arising from their sales/ purchase contract.

POLAR®

Copyright © 2003-2004
Polar Electro Oy, FIN-90440 KEMPELE, Finland

All rights reserved. No part of this manual may be used or reproduced in any form or by any means without prior written permission of Polar Electro Oy.

The names and logos marked with a TM symbol in this user's manual or in the package of this product are trademarks of Polar Electro Oy. The names and logos

marked with a ® symbol in this user's manual or in the package of this product are registered trademarks of Polar Electro Oy.



This CE marking shows compliance of this product with Directive 93/42/EEC.

Polar Electro Oy is a ISO 9001:2000 certified company.

This product is protected by one or several of the following patents:

FI 88223, DE 4215549, FR 92.06120, GB 2257523, HK 113/1996, SG 9591671-4, US 5491474, FI 88972, DE 9219139.8, FR 92.09150, GB 2258587, HK 306/1996, SG 9592117-7, US 5486818, FI 96380, US 5611346, EP 665947, DE 69414362, FI 68734, DE 3439238, GB 2149514, HK 812/1989, US 4625733, FI 100924, FI 100452, US 5840039, FI 4150, DE 20008882.3, US 6477397, FR 0006504, FI 4069, US 6272365, FR 9907823, GB 2339833, DE 29910633, FI 104463, US 6183422, FI 4157, DE 20008883.1, FR 0006778, WO 97/33512, GB 2326240, FI 107776, US 6327486, US 6277080, US 6361502, US 5719825, US 5848027, EP 1055158.

Other patents pending.

15. DISCLAIMER

The material in this manual is for informational purposes only. The products it describes are subject to change without prior notice, due to the manufacturer's continuous development program.

Polar Electro Inc. / Polar Electro Oy makes no representations or warranties with respect to this manual or with respect to the products described herein. Polar Electro Inc. / Polar Electro Oy shall not be liable for any damages, losses, costs or expenses, direct, indirect or incidental, consequential or special, arising out of, or related to the use of this material or the products described herein.

Bullets in the display during exercise	63	Reset the total calorie consumption count	66
Button sound ON/OFF	62	Reset the total exercising time	66
Buttons and their functions	48	Return to Time of day display	49
Care and maintenance	71	See your heart rate	62
Check exercise information	63	See your progress during exercise	63
Check your heart rate limits	63	Setting activity level	51
Disclaimer	80	Setting age	51
Exercise reminder at full hours	49	Setting heart rate limits manually	52
Fitness test settings	67	Setting height	51
Frequently asked questions	75	Setting measurement options	52
Functions during heart rate measurement	62	Setting sex	51
Glossary	77	Setting the alarm	55
Guarantee	79	Setting the date	56
Illuminate the display	62	Setting the time of day	56
Latest OwnIndex	66	Setting user 1 or 2	50
Modes and their functions	49	Setting user information	50
OwnCal calorie consumption	54	Setting weight	51
OwnCal calorie counter ON/OFF	53	Sound signal ON/OFF	62
OwnIndex	67	Start the heart rate measurement	58
OwnZone	59	Start/stop the stopwatch	62
OwnZone ON/OFF	52	Stop the heart rate measurement	64
Polar Fitness Test ON/OFF	54	Subtract the total calorie expenditure amount	65
Precautions	73	Swap the heart rate limits	63
Recalling training information	65	Swap the user	62
Reset stopwatch	60, 63	Technical specifications	78